

1.10024.0001

1.10024.0007

## MQuant® Ammonium Test

 $\text{NH}_4^+$ 

### 1. Method

Ammonium ions react with Neßler's reagent to form a yellow-brown compound. The concentration of ammonium is measured **semiquantitatively** by visual comparison of the reaction zone of the test strip with the fields of a color scale.

### 2. Measuring range and number of determinations

Measuring range / color-scale graduation	Number of determinations
10 - 30 - 60-100-200- <b>400 mg/l <math>\text{NH}_4^+</math></b>	100
8 - 23 - 47 - 78 - 155- <b>310 mg/l <math>\text{NH}_4^+</math>-N</b>	

<sup>1)</sup> for conversion factors see section 8

### 3. Applications

#### Sample material:

Groundwater and surface water

Wastewater

Fertilizers

Process water (e.g. textile industry, plastics industry)

### 4. Influence of foreign substances

This was checked individually in solutions with 200 mg/l  $\text{NH}_4^+$ . The determination is not yet interfered with up to the concentrations of foreign substances given in the table. Cumulative effects were not checked; such effects can, however, not be excluded.

Concentrations of foreign substances in mg/l					
Al <sup>3+</sup>	1000	Cu <sup>2+</sup>	1000	Ni <sup>2+</sup>	100
Ca <sup>2+</sup>	100	Fe <sup>2+</sup>	<b>10</b>	NO <sub>2</sub> <sup>-</sup>	1000
Cl <sup>-</sup>	1000	Fe <sup>3+</sup>	1000	NO <sub>3</sub> <sup>-</sup>	1000
CN <sup>-</sup>	<b>10</b>	K <sup>+</sup>	1000	PO <sub>4</sub> <sup>3-</sup>	1000
Cr <sup>3+</sup>	100	Mg <sup>2+</sup>	1000	S <sub>2</sub> O <sub>3</sub> <sup>2-</sup>	1000
CrO <sub>4</sub> <sup>2-</sup>	1000	Mn <sup>2+</sup>	<b>10</b>		

### 5. Reagents and auxiliaries

#### Please note the warnings on the packaging materials!

The test strips and the test reagent are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

#### Package contents:

Tube containing 100 test strips  
2 bottles of reagent NH<sub>4</sub>-1  
1 test vessel

#### Other reagents:

Ammonium standard solution Certipur®,  
1000 mg/l  $\text{NH}_4^+$ , Cat. No. 119812

### 6. Preparation

Samples containing more than 400 mg/l  $\text{NH}_4^+$  must be diluted with distilled water.

### 7. Procedure

Rinse the test vessel several times with the pre-treated sample.

Pretreated sample (15 - 25 °C)	5 ml	Fill the test vessel to the 5-ml mark.
Reagent NH <sub>4</sub> -1	10 drops <sup>1)</sup>	Add and swirl.

Immerse the reaction zone of the test strip in the measurement sample **for 3 sec.**

Allow excess liquid to run off via the long edge of the strip onto an absorbent paper towel and **after 10 sec** determine with which color field on the label the color of the reaction zone coincides most exactly.

Read off the corresponding result in mg/l  $\text{NH}_4^+$  or NH<sub>4</sub>-N.

<sup>1)</sup> Hold the bottle vertically while adding the reagent!

#### Notes on the measurement:

- The color of the reaction zone may continue to change after the specified reaction time has elapsed. This must not be considered in the measurement.
- If the color of the reaction zone is equal to or more intense than the darkest color on the scale, repeat the measurement using **fresh**, diluted samples until a value of less than 400 mg/l  $\text{NH}_4^+$  is obtained.

Concerning the result of the analysis, the dilution (see also section 6) must be taken into account:

Result of analysis = measurement value x dilution factor

### 8. Conversions

Units required	= units given x conversion factor
mg/l NH <sub>4</sub> -N	mg/l $\text{NH}_4^+$ 0.776
mg/l $\text{NH}_4^+$	mg/l NH <sub>4</sub> -N 1.29

### 9. Method control

To check test strips, test reagent, and handling: Dilute the ammonium standard solution with distilled water to 100 mg/l  $\text{NH}_4^+$  and analyze as described in section 7.

Additional notes see under [www.qa-test-kits.com](http://www.qa-test-kits.com).

### 10. Notes

- Reclose** the reagent bottle and **the tube containing the test strips immediately after use.**
- Rinse the test vessel **with distilled water only.**

