FIELD SCOUT® SOILSTIKTM PH METER OPERATION MANUAL

CATALOG #2105



Spectrum ® Technologies, Inc.

CONTENTS

General Overview3
LCD Display4
LCD Display Messages5
Meter Components6
Calibration7
Changing the Temperature Units8
Taking Liquid/Soil Measurements9
Storing/Recalling Measurements10
Sensor Replacement11
Battery Replacement12
Meter Handling/Troubleshooting13
Specifications14
Warranty15

This manual will familiarize you with the features and operation of your new meter. Please read this manual thoroughly before using your meter. For customer support or to place an order, contact Spectrum Technologies, Inc between 7:30 a.m. and 5:30 p.m. CST.

Phone: (800) 248-8873 or (815) 436-4440 Fax: (815) 436-4460

E-mail: info@specmeters.com **Website:** www.specmeters.com

GENERAL OVERVIEW

Congratulations on the purchase of your Field Scout[™] SoilStik[™] Pro Meter. This manual describes how to use your pH meter and keep it working accurately. Read the manual thoroughly in order to make effective use of your meter.


The SoilStik delivers high quality answers, with an accuracy of +/- 0.01 pH units. This self-contained digital meter allows you to test the pH levels in water, soil, and other liquids. The replaceable sensor makes the measurement of small samples much more convenient. The sensor provides a visual indication of when replacement is due.

The meter has a two-point, automatic calibration (pH 4 and 7), with a range of pH 2.0-12.0. The display will show your results to a resolution of 0.01 pH units.

Getting Started

- New meters are shipped with a battery insulating strip. Remove this strip before attempting to power up the meter. The battery compartment is accessed by removing the cap at the top of the meter.
- Remove the cap from the bottom of the SoilStik to expose the electrode glass surface and reference junction.
- Before first use, or after extended storage, soak the electrode (with its cap removed) in pH 4 solution for about 10 minutes.

LCD DISPLAY



- 1. Bar graph reading
- 2. Measurement reading
- 3. BAT (low battery) and HOLD (data hold) indicators
- 4. Temperature display
- 5. Bar graph scale designations
- 6. Units of measure
- 7. Calibration indicators
- 8. RENEW and CAL indicators

LCD DISPLAY MESSAGES

In addition to displaying the pH and temperature readings, the LCD of the SoilStik can also display:

- CAL reminder

If the SoilStik is powered up in pH mode 15 times without being calibrated, the **CAL** reminder symbol appears on the LCD indicating that the SoilStik may require calibration. Some applications may require recalibration of the electrode more frequently than others. The **CAL** display is simply a reminder and will turn off when the pH electrode is recalibrated.

- **RENEW** reminder

The **RENEW** symbol appears when the output of the pH electrode fails an internal diagnostic test. A flashing **RENEW** warning indicates that the probe <u>may</u> be nearing the end of its useful life. Try to clean and recalibrate the meter. If the RENEW symbol remains flashing, Remove the battery holder, remove and replace the batteries and re-insert the battery holder (see p. 12). If this does not cause the **RENEW** symbol to disappear, replace the electrode.


- HOLD

While in measurement mode, pressing the **MODE** button causes the meter to hold the value on the LCD. The reading is also stored in short-term memory (see p. 10)

- BAT

Indicates that battery power is low. To preserve battery life, the meter automatically powers off after 10 minutes of inactivity.

METER COMPONENTS



- 1. Battery compartment cap
- 2. LCD Display
- 3. MODE/HOLD button
- 4. CAL/RECALL button
- 5. **ON/OFF** button
- 6. Electrode collar
- 7. pH electrode

(Electrode cap is not shown)

CALIBRATION

The SoilStik automatically recognizes buffer solution pH's of 4, 7, and 10. Buffer solutions of pH 4 and 7 are included with the meter. Although it is possible to calibrate to all 3 values, a two-point calibration with a buffer of 7 plus 4 or 10 (whichever is nearest to the expected sample value) is always recommended. A one-point calibration (with a value closest to the expected sample value) is also valid. For best accuracy, always calibrate at the sample temperature.

Calibration Procedure

- 1. Put a small amount of buffer solution in a small plastic cup (provided). For multi-point calibrations, pH 7 should be calibrated first, followed by 4 and/or 10. Place the electrode into the buffer solution and momentarily press the **CAL/RECALL** button.
- 2. The SoilStik will automatically recognize the calibration solution. After the **CAL/RECALL** button is released, the LCD will display **CAL/RECALL**. The pH calibration value will then flash on the LCD.
- 3. When calibration is complete, the SoilStik automatically displays the **END** symbol and returns to normal operation mode. Note: If the calibration solution is more than 1 pH unit different from pH 4, 7, or 10, the SoilStik will assume an error has been made, abort the calibration and return to normal operation mode without changing any internal calibration constants.
- 4. The appropriate circled indicator (, , or) will appear on the LCD when a calibration has been completed with a given calibration buffer solution. The calibration constants are stored until a new calibration is performed.
- 5. For a two- or three point-calibration, repeat steps 1 to 4. As additional calibrations are performed, additional circled indicators will appear on the LCD.

When the meter is powered off, the calibration constants are retained. Therefore, the circled indicators for all completed calibrations will remain illuminated. However, when a new calibration is done, only the circled indicator for that pH level will be displayed. If a multiple-point calibration is desired, the calibration must be performed for the other pH levels as well.

Notes:

- If the batteries are removed or go dead, any usercreated calibration data will be lost and the meter will return to the factory defaults.
- Always turn the meter off and then on before calibrating to allow sufficient time to complete the calibrations during one power cycle, especially for multi-point calibrations. If the meter powers off automatically during a calibration, any completed calibrations will be saved. But, after the meter is powered up again, any new multi-point calibration will have to be done in entirety.

CHANGING THE TEMPERATURE UNITS

Press and hold the **CAL/RECALL** button for approximately 3 seconds. The °C or °F symbol will change first and the numerical temperature value will change *after* the button is released. If the meter accidentally transitions into **CAL** mode, the **CAL** symbol will appear on the LCD. Simply turn the SoilStik off and start again.

TAKING LIQUID AND SOIL MEASUREMENTS

Liquids Readings

When the electrode is placed in a solution, the main display and bar graph indicate the pH reading while the lower display reads temperature (readings flash until they have stabilized).

The curved bar graph is centered at pH 7. At a pH of 7, no bars will be visible on the bar graph. When the pH > 7, the bars will extend to the right. When the pH < 7, the bars will extend to the to the left.

Soil Slurry Readings

- 1. Add 1 part soil to 1 part distilled water.
- 2. Mix for 30 seconds.
- 3. Allow 1 to 2 minutes for soil slurry to come into equilibrium
- 4. Submerge the sensor into the soil slurry and record the result.
- 5. Wash the sensor with distilled water.

Direct Soil Readings

The SoilStik can also be used to take measurements on soil cores. Moisten the surface of the core with distilled water and press the electrode onto the moist soil. Be sure that the clear sensor and reference electrode are both in contact with the soil. Rinse and clean the sensor well before taking additional readings.

STORING AND RECALLING MEASUREMENTS

Storing readings in short-term memory

While in measurement mode, pressing the MODE/HOLD button will freeze the value on the LCD and store the number in memory. The LCD will briefly display the memory location number and then the value stored. The HOLD symbol will appear above the temperature reading. Momentarily press MODE/HOLD again to return to normal operation.

The short-term memory can hold 15 readings. After 15 readings have been saved, the SoilStik will return to memory location 1 and start overwriting existing data with newly stored data.

Recalling stored readings

Check that the **HOLD** symbol is not displayed. If it is, exit the **HOLD** function by momentarily pressing the **MODE/HOLD** button.

To view stored readings, momentarily press the CAL/RECALL button and then press the MODE/HOLD button while the CAL symbol is displayed. The storage location number (1 through 15) will flash. If the meter accidentally transitions into CAL mode (causing the display to flash), press the CAL/RECALL button again to exit and start over.

The last reading stored will be displayed first. To advance through the stored readings, momentarily press the MODE/HOLD button. The location number is displayed first, followed by the reading stored in that location. To exit recall mode, momentarily press the CAL/RECALL button and the SoilStik will return to normal operation.

Note: If the batteries are removed or go dead, any stored readings will be lost.

SENSOR REPLACEMENT

The SoilStik is shipped with an electrode attached. Electrode life is limited and is dependent on (among other factors) frequency of use and care. If the electrode needs to be replaced, follow these steps for removing and connecting electrodes.


- 1. To remove an electrode, unscrew and completely remove the electrode retaining collar.
- 2. Gently rock the electrode from side to side, pulling it away from the meter, until it disconnects.
- 3. To attach an electrode, carefully plug the electrode into the meter socket (note that the electrode connector is keyed, ensuring proper connection).
- 4. Secure the electrode in place by tightly turning the collar in place. (a rubber gasket seals the electrode with the meter).

BATTERY REPLACEMENT

The batteries need to be replaced when the LCD shows **BAT** symbol while the power is turned on or if the LCD shows nothing when you press the **ON/OFF** button.

Battery Replacement Procedure:

- 1. Twist off the battery compartment cap.
- 2. Slide the black battery holder out of the meter chassis.
- 3. Replace the four (4) CR-2032 batteries. Take care that the positive side is facing up.
- 4. Return the battery holder. Be sure to align the 2 metal strips.
- 5. Replace the battery compartment cap.



METER HANDLING/ TROUBLESHOOTING

- Be sure the sponge inside the electrode cap is hydrated at all times with pH buffer buffer solutions.
- White KCI crystals may be present in the cap. These crystals will dissolve while the sensor is soaked or they can be simply rinsed off with tap water
- If the unit appears to be locked (display frozen) it is possible that the meter is in **Hold** mode (see p. 12) Simply press the **MODE/HOLD** button again or turn the meter off and restart if the display appears frozen. If this fails to revive the, meter, remove the batteries, push the ON button for 3 seconds and then reinsert the batteries.
- If the meter has difficulty calibrating or maintaining a stable reading, the sensor may need to be replaced.
- Accuracy will be improved if the standard solution is at the same temperature as that of sample being measured.

SPECIFICATIONS

Calibration: Two-point, pH 4 and 7

Display: Digital LCD with bar graph.

Range: pH 0.00 - pH 12.00

Accuracy: +/- 0.01 pH

Resolution: +/- 0.01 pH

Operating 32 to $122^{\circ}F$ (0 to $50^{\circ}C$) / < 80% RH **Temperature:**

•

Temperature 23 to 194°F (-5 to 90°C) **Range**:

Temperature 0.1° up to 99.9 then 1° thereafter **Resolution:**

Temperature $\pm 1.8^{\circ}\text{F} / 1^{\circ}\text{C}$ [from 23 to 122°F **Accuracy:** (-5 to 50°C)]

± 5.4°F / 3°C [from 122 to 194°F

(50 to 90°C)]

Short-term 15 tagged (numbered) readings memory:

Battery: Four (4)CR2032 button batteries

Auto Power-off: After 10 minutes of inactivity

Dimensions: 35.6 x 172.7 x 40.6 mm

(1.4 x 6.8 x 1.6 inches)

Weight: 110g (3.85oz)



WARRANTY

This product has been brought to you having passed severe quality control and inspections. Should any trouble occur during the course of normal use, the meter shall be repaired or replaced free of charge in accordance with the stipulations laid down herein. The term of this warranty shall be for one year from date of purchase. This warranty excludes batteries, sensor and accessories.

Warranty Stipulations:

- 1. The product shall be repaired or replaced free of charge should any trouble occur during the course of normal use if returned within the warranty period (one year from date of purchase). In which event, contact the dealer of purchase. Return the meter with proof of date of purchase.
- 2. Expenses shall be incurred in the following instances within the warranty period. (Costs such as postage shall be born by the customer)
- a) When the date of purchase and store name is not written on the warranty.
- b) When trouble or damage has been incurred due to misuse, abuse, and/or improper handling.
- c) When the meter has been repaired, modified and dismantled by persons other than the designated agent or service shop.
- d) In the event of changes in external appearance such as scratches or dirt caused during use or battery fluid leakage.
- e) In the event of unsuitable movement, dropping or accidents such as fire, earthquakes, floods or a burglary.
- f) When replacing consumables and accessories.
- g) When cause of trouble lies not in the meter itself.
- h) When this warranty is not shown and when necessary particulars have not been written in the warranty.

Our obligation under this warranty is to repair or replace the meter free of charge in accordance with the conditions laid down herein. Accordingly, this warranty doesn't limit your specific legal rights.

Product Return

If for any reason you are not satisfied, or the meter has failed and you need to return the product for service, you need to contact Spectrum Technologies, Inc.

Before returning a failed unit, you must obtain a Returned Goods Authorization (RGA) number from Spectrum Technologies. You must ship the product(s), properly packaged against further damage, back to Spectrum Technologies (at your expense) with the RGA number marked clearly on the outside of the package. Spectrum Technologies is not responsible for any package that is returned without a valid RGA number, or for the loss of the package by any shipping company.

Service and Support

For technical service and support call your distributor or Spectrum Technologies, Inc.

When calling for technical support, have a detailed explanation of the problem that you are experiencing. The more information you can provide, the faster and easier a technical support person will be able to assist you.

Spectrum Technologies, Inc. 12360 S Industrial Dr East Plainfield, IL 60585 (800) 248-8873 or (815) 436-4440 Fax (815) 436-4460

E-Mail: info@specmeters.com Website: www.specmeters.com